

## Tbx3 Rabbit mAb

Catalog No: #56724

Package Size: #56724-1 50ul #56724-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Tbx3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP FC
Species Reactivity	Human Mouse Rat
Specificity	Tbx3 Antibody detects endogenous levels of total Tbx3
Immunogen Description	A synthesized peptide derived from human Tbx3
Other Names	T box 3; T-box protein 3; Tbx3; TBX3 ISO; UMS; XHL;
Accession No.	Uniprot:O15119
Uniprot	O15119
Calculated MW	85kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IP:1:50FC:1:100

## Images



Western blot analysis of Tbx3 expression in MCF7 cell lysate.

## Product Description

Transcriptional repressor involved in developmental processes. Probably plays a role in limb pattern formation.

---

Note: This product is for in vitro research use only