

## Thromboxane synthase Rabbit mAb

Catalog No: #56728

Package Size: #56728-1 50ul #56728-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

|                       |  |
|-----------------------|--|
| Product Name          | Thromboxane synthase Rabbit mAb  |
| Host Species          | Rabbit   |
| Clonality             | Monoclonal   |
| Isotype               | Rabbit IgG   |
| Purification          | Affinity-chromatography  |
| Applications          | WB   |
| Species Reactivity    | Human Mouse  |
| Specificity           | Thromboxane synthase Antibody detects endogenous levels of total Thromboxane synthase              |
| Immunogen Description | A synthesized peptide derived from human Thromboxane synthase                                      |
| Other Names           | TS; TXS; CYP5; THAS; TXAS; CYP5A1; GHOSAL; BDPLT14;  |
| Accession No.         | Uniprot:P24557   |
| Uniprot               | P24557   |
| Calculated MW         | 54kDa  |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.                   |

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of Thromboxane synthase expression in Human fetal liver lysate.

## Product Description

This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids.

---

Note: This product is for in vitro research use only