

GLUD1 Rabbit mAb

Catalog No: #56733

Package Size: #56733-1 50ul #56733-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GLUD1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	GLUD1 Antibody detects endogenous levels of total GLUD1
Immunogen Description	A synthesized peptide derived from human GLUD1
Other Names	GDH; Gdh-X; GDH1; GLUD; Glud1; Glud1a; Glud1b; Gludl; MRG-2;
Accession No.	Uniprot:P00367
Uniprot	P00367
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:500

Images



Western blot analysis of GLUD1 expression in HepG2 cell lysate.

Product Description

May be involved in learning and memory reactions by increasing the turnover of the excitatory neurotransmitter glutamate.

Note: This product is for in vitro research use only