

LOX12 Rabbit mAb

Catalog No: #56736

Package Size: #56736-1 50ul #56736-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LOX12 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Rat
Specificity	LOX12 Antibody detects endogenous levels of total LOX12
Immunogen Description	A synthesized peptide derived from human LOX12
Other Names	12 LOX; 12LO; 12S LOX; 12S-lipoxygenase; 12S-LOX; Alox12; LOG12; P-12LO;
Accession No.	Uniprot:P18054
Uniprot	P18054
Calculated MW	70kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of LOX12 expression in A431 cell lysate.

Product Description

Oxygenase and 14,15-leukotriene A4 synthase activity.

Note: This product is for in vitro research use only