

HOXB4 Rabbit mAb

Catalog No: #56739

Package Size: #56739-1 50ul #56739-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HOXB4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	HOXB4 Antibody detects endogenous levels of total HOXB4
Immunogen Description	A synthesized peptide derived from human HOXB4
Other Names	HOX2; Hox2.6; Hox2F; hoxb4;
Accession No.	Uniprot:P17483
Uniprot	P17483
Calculated MW	34kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of HOXB4 expression in Jurkat cell lysate.

Product Description

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Note: This product is for in vitro research use only