

DNAJA1 Rabbit mAb

Catalog No: #56764

Package Size: #56764-1 50ul #56764-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | DNAJA1 Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IP FC |
| Species Reactivity | Human Mouse Rat |
| Specificity | DNAJA1 Antibody detects endogenous levels of total DNAJA1 |
| Immunogen Description | A synthesized peptide derived from human DNAJA1 |
| Other Names | DJ2; DjA1; DNAJ2; Dnaja1; HDJ2; HSDJ; HSJ2; HSPF4; NEDD7; |
| Accession No. | Uniprot:P31689 |
| Uniprot | P31689 |
| Calculated MW | 45kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:1000~1:5000IP:1:50FC:1:50

Images



Western blot analysis of DNAJA1 expression in Jurkat cell lysate.

Product Description

Co-chaperone of Hsc70. Seems to play a role in protein import into mitochondria.

Note: This product is for in vitro research use only