

DCTN3 Rabbit mAb

Catalog No: #56772

Package Size: #56772-1 50ul #56772-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DCTN3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	DCTN3 Antibody detects endogenous levels of total DCTN3
Immunogen Description	A synthesized peptide derived from human DCTN3
Other Names	DCNT22; DCTN22; Dctn3; dynactin 3 (p22); dynactin light chain; Dynactin subunit 3; p22;
Accession No.	Uniprot:O75935
Uniprot	O75935
Calculated MW	21kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of DCTN3 expression in SH-SY5Y cell lysate.

Product Description

Together with dynein may be involved in spindle assembly and cytokinesis.

Note: This product is for in vitro research use only