

DP1 Rabbit mAb

Catalog No: #56773

Package Size: #56773-1 50ul #56773-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	DP1 Antibody detects endogenous levels of total DP1
Immunogen Description	A synthesized peptide derived from human DP1
Other Names	DP1; DRTF1; TFDP1;
Accession No.	Uniprot:Q14186
Uniprot	Q14186
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000FC:1:50

Images



Western blot analysis of DP1 expression in Raji cell lysate.

Product Description

Can stimulate E2F-dependent transcription. Binds DNA cooperatively with E2F family members through the E2 recognition site, 5'-TTTC[CG]CGC-3', found in the promoter region of a number of genes whose products are involved in cell cycle regulation or in DNA replication.

Note: This product is for in vitro research use only