

TMOD2 Rabbit mAb

Catalog No: #56777

Package Size: #56777-1 50ul #56777-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TMOD2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	TMOD2 Antibody detects endogenous levels of total TMOD2
Immunogen Description	A synthesized peptide derived from human TMOD2
Other Names	Tropomodulin 2; NTMOD; TMOD2;
Accession No.	Uniprot:Q9NZR1
Uniprot	Q9NZR1
Calculated MW	40kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of TMOD2 expression in SH-SY5Y cell lysate.

Product Description

Blocks the elongation and depolymerization of the actin filaments at the pointed end. The Tmod/TM complex contributes to the formation of the short actin protofilament, which in turn defines the geometry of the membrane skeleton (By similarity).

Note: This product is for in vitro research use only