

ANGPTL3 Rabbit mAb

Catalog No: #56782

Package Size: #56782-1 50ul #56782-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ANGPTL3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	ANGPTL3 Antibody detects endogenous levels of total ANGPTL3
Immunogen Description	A synthesized peptide derived from human ANGPTL3
Other Names	ANG5; Angiotensin II type 1 receptor 1; Angiotensin II type 1 receptor 2; ANGPT5; ANGPTL3; ANL3; FHBL2;
Accession No.	Uniprot:Q9Y5C1
Uniprot	Q9Y5C1
Calculated MW	54kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 ICC/IF: 1:50~1:200

Images



Western blot analysis of ANGPTL3 expression in A375 cell lysate.

Product Description

Acts in part as a hepatokine that is involved in regulation of lipid and glucose metabolism. Proposed to play a role in the trafficking of energy substrates to either storage or oxidative tissues in response to food intake (By similarity). Has a stimulatory effect on plasma triglycerides (TG), which is achieved by suppressing plasma TG clearance via inhibition of LPL activity.

Note: This product is for in vitro research use only