

Desmoglein 3 Rabbit mAb

Catalog No: #56784

Package Size: #56784-1 50ul #56784-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Desmoglein 3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Desmoglein 3 Antibody detects endogenous levels of total Desmoglein 3
Immunogen Description	A synthesized peptide derived from human Desmoglein 3
Other Names	CDHF6; Desmoglein3; DSG3; PVA;
Accession No.	Uniprot:P32926
Uniprot	P32926
Calculated MW	145kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of Desmoglein 3 expression in A431 cell lysate.

Product Description

Component of intercellular desmosome junctions. Involved in the interaction of plaque proteins and intermediate filaments mediating cell-cell adhesion.

Note: This product is for in vitro research use only