

FABPI Rabbit mAb

Catalog No: #56789

Package Size: #56789-1 50ul #56789-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FABPI Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	FABPI Antibody detects endogenous levels of total FABPI
Immunogen Description	A synthesized peptide derived from human FABPI
Other Names	FABP2; FABPI; IFABP;
Accession No.	Uniprot:P12104
Uniprot	P12104
Calculated MW	15kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of intestinal FABP expression in Human small intestine lysate.

Product Description

FABP are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters. FABP2 is probably involved in triglyceride-rich lipoprotein synthesis.

Note: This product is for in vitro research use only