

## Inhibin beta B Rabbit mAb

Catalog No: #56790

Package Size: #56790-1 50ul #56790-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Inhibin beta B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Inhibin beta B Antibody detects endogenous levels of total Inhibin beta B
Immunogen Description	A synthesized peptide derived from human Inhibin beta B
Other Names	Activin beta B; INHBB;
Accession No.	Uniprot:P09529
Uniprot	P09529
Calculated MW	52kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of Inhibin beta B expression in Mouse brain lysate.

## Product Description

Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary gland. Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival, embryonic axial development or bone growth, depending on their subunit composition.

---

Note: This product is for in vitro research use only