

NAIP Rabbit mAb

Catalog No: #56792

Package Size: #56792-1 50ul #56792-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NAIP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	NAIP Antibody detects endogenous levels of total NAIP
Immunogen Description	A synthesized peptide derived from human NAIP
Other Names	BIRC1; Birc1a; NAIP; Naip1; NLRB1; psiNAIP;
Accession No.	Uniprot:Q13075
Uniprot	Q13075
Calculated MW	160kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of NAIP expression in Human fetal kidney lysate.

Product Description

Prevents motor-neuron apoptosis induced by a variety of signals. Possible role in the prevention of spinal muscular atrophy that seems to be caused by inappropriate persistence of motor-neuron apoptosis: mutated or deleted forms of NAIP have been found in individuals with severe spinal muscular atrophy.

Note: This product is for in vitro research use only