

SSTR5 Rabbit mAb

Catalog No: #56796

Package Size: #56796-1 50ul #56796-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SSTR5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	SSTR5 Antibody detects endogenous levels of total SSTR5
Immunogen Description	A synthesized peptide derived from human SSTR5
Other Names	SSTR5; SS5R; Sstr5;
Accession No.	Uniprot:P35346
Uniprot	P35346
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000FC:1:100

Images



Western blot analysis of Somatostatin Receptor 5 expression in HeLa cell lysate.

Product Description

Receptor for somatostatin 28 and to a lesser extent for somatostatin-14. The activity of this receptor is mediated by G proteins which inhibit adenylyl cyclase.

Note: This product is for in vitro research use only