

MFF Rabbit mAb

Catalog No: #56797

Package Size: #56797-1 50ul #56797-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MFF Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	MFF Antibody detects endogenous levels of total MFF
Immunogen Description	A synthesized peptide derived from human MFF
Other Names	C2orf33; GL004; Mff;
Accession No.	Uniprot:Q9GZY8
Uniprot	Q9GZY8
Calculated MW	34kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of MFF expression in Human fetal muscle lysate.

Product Description

Plays a role in mitochondrial and peroxisomal fission.

Note: This product is for in vitro research use only