

## ESRRG Rabbit mAb

Catalog No: #56799

Package Size: #56799-1 50ul #56799-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ESRRG Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	ESRRG Antibody detects endogenous levels of total ESRRG
Immunogen Description	A synthesized peptide derived from human ESRRG
Other Names	ERR 3; ERR gamma 2; ERR3; ERRG2; ERRgamma; ESRRG; NR3B3;
Accession No.	Uniprot:P62508
Uniprot	P62508
Calculated MW	51kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of Estrogen Related Receptor gamma expression in Human fetal heart lysate.

## Product Description

Orphan receptor that acts as transcription activator in the absence of bound ligand. Binds specifically to an estrogen response element and activates reporter genes controlled by estrogen response elements.

---

Note: This product is for in vitro research use only