

## TFPI Rabbit mAb

Catalog No: #56802

Package Size: #56802-1 50ul #56802-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TFPI Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	TFPI Antibody detects endogenous levels of total TFPI
Immunogen Description	A synthesized peptide derived from human TFPI
Other Names	EPI; LACI; TFI; TFPI1;
Accession No.	Uniprot:P10646
Uniprot	P10646
Calculated MW	35-45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of TFPI expression in U87-MG cell lysate.

## Product Description

Inhibits factor X (X(a)) directly and, in a Xa-dependent way, inhibits VIIa/tissue factor activity, presumably by forming a quaternary Xa/LACI/VIIa/TF complex. It possesses an antithrombotic action and also the ability to associate with lipoproteins in plasma.

---

Note: This product is for in vitro research use only