

## Oncostatin M Rabbit mAb

Catalog No: #56803

Package Size: #56803-1 50ul #56803-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Oncostatin M Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Oncostatin M Antibody detects endogenous levels of total Oncostatin M
Immunogen Description	A synthesized peptide derived from human Oncostatin M
Other Names	Oncostatin M; OSM;
Accession No.	Uniprot:P13725
Uniprot	P13725
Calculated MW	22kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of Oncostatin M expression in Jurkat cell lysate.

## Product Description

Growth regulator. Inhibits the proliferation of a number of tumor cell lines. Stimulates proliferation of AIDS-KS cells. It regulates cytokine production, including IL-6, G-CSF and GM-CSF from endothelial cells.

---

Note: This product is for in vitro research use only