

## Wnt5b Rabbit mAb

Catalog No: #56808

Package Size: #56808-1 50ul #56808-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Wnt5b Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	Wnt5b Antibody detects endogenous levels of total Wnt5b
Immunogen Description	A synthesized peptide derived from human Wnt5b
Other Names	Wnt5b;
Accession No.	Uniprot:Q9H1J7
Uniprot	Q9H1J7
Calculated MW	45kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000FC:1:50

## Images



Western blot analysis of Wnt5b expression in HepG2 cell lysate.

## Product Description

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.

---

Note: This product is for in vitro research use only