

## HLA F Rabbit mAb

Catalog No: #56816

Package Size: #56816-1 50ul #56816-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	HLA F Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	HLA F Antibody detects endogenous levels of total HLA F
Immunogen Description	A synthesized peptide derived from human HLA F
Other Names	HLAF;
Accession No.	Uniprot:P30511
Uniprot	P30511
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200

## Images



Western blot analysis of HLA F expression in JAR cell lysate.

## Product Description

Involved in the presentation of foreign antigens to the immune system.

---

Note: This product is for in vitro research use only