

Neurocan Rabbit mAb

Catalog No: #56817

Package Size: #56817-1 50ul #56817-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Neurocan Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	Neurocan Antibody detects endogenous levels of total Neurocan
Immunogen Description	A synthesized peptide derived from human Neurocan
Other Names	Cspg3; Ncan; NEUR; Neurocan core protein;
Accession No.	Uniprot:O14594
Uniprot	O14594
Calculated MW	130,270kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of Neurocan expression in Fetal brain cell lysate.

Product Description

May modulate neuronal adhesion and neurite growth during development by binding to neural cell adhesion molecules (NG-CAM and N-CAM).
Chondroitin sulfate proteoglycan; binds to hyaluronic acid.

Note: This product is for in vitro research use only