

LMO4 Rabbit mAb

Catalog No: #56820

Package Size: #56820-1 50ul #56820-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	LMO4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse
Specificity	LMO4 Antibody detects endogenous levels of total LMO4
Immunogen Description	A synthesized peptide derived from human LMO4
Other Names	LIM domain only protein 4; LMO4;;
Accession No.	Uniprot:P61968
Uniprot	P61968
Calculated MW	17kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IP: 1:50

Images



Western blot analysis of LMO4 expression in Jurkat cell lysate.

Product Description

Probable transcriptional factor.

Note: This product is for in vitro research use only