

PDLIM1 Rabbit mAb

Catalog No: #56823

Package Size: #56823-1 50ul #56823-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PDLIM1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	PDLIM1 Antibody detects endogenous levels of total PDLIM1
Immunogen Description	A synthesized peptide derived from human PDLIM1
Other Names	CLIM1; CLP36; Elfin; Enigma; hCLIM1; Pdlim1;
Accession No.	Uniprot:O00151
Uniprot	O00151
Calculated MW	36kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000 ICC/IF:1:50~1:200 FC:1:50

Images



Western blot analysis of PDLIM1 expression in 293T cell lysate.

Product Description

Cytoskeletal protein that may act as an adapter that brings other proteins (like kinases) to the cytoskeleton.

Note: This product is for in vitro research use only