

## Eph receptor B3 Rabbit mAb

Catalog No: #56827

Package Size: #56827-1 50ul #56827-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Eph receptor B3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	Eph receptor B3 Antibody detects endogenous levels of total Eph receptor B3
Immunogen Description	A synthesized peptide derived from human Eph receptor B3
Other Names	Cek10; EK2; Embryonic kinase 2; ephb3; ETK2; hEK2; Mdk5; Sek4; TYRO6;
Accession No.	Uniprot:P54753
Uniprot	P54753
Calculated MW	110kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:50FC:1:50

## Images



Western blot analysis of Eph receptor B3 expression in MOLT4 cell lysate.

## Product Description

Receptor for members of the ephrin-B family. Binds to ephrin-B1 and -B2.

---

Note: This product is for in vitro research use only