

DPYD Rabbit mAb

Catalog No: #56829

Package Size: #56829-1 50ul #56829-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DPYD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human
Specificity	DPYD Antibody detects endogenous levels of total DPYD
Immunogen Description	A synthesized peptide derived from human DPYD
Other Names	DHP; DHPDHase; DPD; DPYD;
Accession No.	Uniprot:Q12882
Uniprot	Q12882
Calculated MW	110kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000 IHC: 1:50~1:200 ICC/IF: 1:50~1:200

Images



Western blot analysis of DPYD expression in HeLa cell lysate.

Product Description

Involved in pyrimidine base degradation. Catalyzes the reduction of uracil and thymine. Also involved the degradation of the chemotherapeutic drug 5-fluorouracil.

Note: This product is for in vitro research use only