

GAA Rabbit mAb

Catalog No: #56835

Package Size: #56835-1 50ul #56835-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GAA Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	GAA Antibody detects endogenous levels of total GAA
Immunogen Description	A synthesized peptide derived from human GAA
Other Names	Acid maltase; GAA; Glucosidase alpha acid; LYAG;
Accession No.	Uniprot:P10253
Uniprot	P10253
Calculated MW	120,100,75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of GAA expression in Human fetal liver lysate.

Product Description

Essential for the degradation of glycogen to glucose in lysosomes.

Note: This product is for in vitro research use only