

## NFYA Rabbit mAb

Catalog No: #56839

Package Size: #56839-1 50ul #56839-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	NFYA Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	NFYA Antibody detects endogenous levels of total NFYA
Immunogen Description	A synthesized peptide derived from human NFYA
Other Names	CBFA; CBFB; HAP2; NFYA; Sez10;
Accession No.	Uniprot:P23511
Uniprot	P23511
Calculated MW	37kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IC/IF:1:50~1:200IP:1:50FC:1:50

## Images



Western blot analysis of NFYA expression in 293T cell lysate.

## Product Description

Stimulates the transcription of various genes by recognizing and binding to a CCAAT motif in promoters, for example in type 1 collagen, albumin and beta-actin genes.

---

Note: This product is for in vitro research use only