

GUSB Rabbit mAb

Catalog No: #56848

Package Size: #56848-1 50ul #56848-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | GUSB Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB IHC ICC/IF FC |
| Species Reactivity | Human |
| Specificity | GUSB Antibody detects endogenous levels of total GUSB |
| Immunogen Description | A synthesized peptide derived from human GUSB |
| Other Names | Beta-glucuronidase; Beta-G1; GUSB; |
| Accession No. | Uniprot:P08236 |
| Uniprot | P08236 |
| Calculated MW | 78kDa |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle. |

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of GUSB expression in K562 cell lysate.

Product Description

Plays an important role in the degradation of dermatan and keratan sulfates.

Note: This product is for in vitro research use only