

## DCTN5 Rabbit mAb

Catalog No: #56852

Package Size: #56852-1 50ul #56852-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	DCTN5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	DCTN5 Antibody detects endogenous levels of total DCTN5
Immunogen Description	A synthesized peptide derived from human DCTN5
Other Names	DCTN5; Dynactin 4; Dynactin 5 (p25);
Accession No.	Uniprot:Q9BTE1
Uniprot	Q9BTE1
Calculated MW	22kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

## Images



Western blot analysis of DCTN5 expression in MCF7 cell lysate.

## Product Description

Member of the pointed-end complex of the dynactin shoulder complex which contains DCTN4, DCTN5 and DCTN6 subunits and ACTR10. Within the complex DCTN6 forms a heterodimer with DCTN5.

---

Note: This product is for in vitro research use only