

IL32 Rabbit mAb

Catalog No: #56854

Package Size: #56854-1 50ul #56854-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IL32 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	IL32 Antibody detects endogenous levels of total IL32
Immunogen Description	A synthesized peptide derived from human IL32
Other Names	IL 32alpha; IL 32beta; IL 32delta; IL 32gamma; IL32; NK4; TAIF; TAIFa; TAIFb; TAIFc; TAIFd;
Accession No.	Uniprot:P24001
Uniprot	P24001
Calculated MW	23kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of IL32 expression in Human tonsil lysate.

Product Description

Cytokine that may play a role in innate and adaptive immune responses. It induces various cytokines such as TNFA/TNF-alpha and IL8. It activates typical cytokine signal pathways of NF-kappa-B and p38 MAPK.

Note: This product is for in vitro research use only