

ETV6 Rabbit mAb

Catalog No: #56856

Package Size: #56856-1 50ul #56856-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ETV6 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	ETV6 Antibody detects endogenous levels of total ETV6
Immunogen Description	A synthesized peptide derived from human ETV6
Other Names	ETV6; Tel; TEL oncogene; TEL/ABL; TEL1;
Accession No.	Uniprot:P41212
Uniprot	P41212
Calculated MW	53kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of ETV6 expression in A431 cell lysate.

Product Description

Transcriptional repressor; binds to the DNA sequence 5'-CCGGAAGT-3'.

Note: This product is for in vitro research use only