

HOXA9 Rabbit mAb

Catalog No: #56861

Package Size: #56861-1 50ul #56861-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HOXA9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	HOXA9 Antibody detects endogenous levels of total HOXA9
Immunogen Description	A synthesized peptide derived from human HOXA9
Other Names	HOX1; ABD-B; HOX1G; HOX1.7;
Accession No.	Uniprot:P31269
Uniprot	P31269
Calculated MW	35kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of HOXA9 expression in 293T cell lysate.

Product Description

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Note: This product is for in vitro research use only