

Dynamin 3 Rabbit mAb

Catalog No: #56865

Package Size: #56865-1 50ul #56865-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Dynamin 3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse
Specificity	Dynamin 3 Antibody detects endogenous levels of total Dynamin 3
Immunogen Description	A synthesized peptide derived from human Dynamin 3
Other Names	DNM3; Dyn3; DynaIII; dynamin 2; Dynamin; Dynamin3; T dynamin;
Accession No.	Uniprot:Q9UQ16
Uniprot	Q9UQ16
Calculated MW	97kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of Dynamin 3 expression in Human fetal brain lysate.

Product Description

Microtubule-associated force-producing protein involved in producing microtubule bundles and able to bind and hydrolyze GTP. Most probably involved in vesicular trafficking processes, in particular endocytosis.

Note: This product is for in vitro research use only