

Complement factor B Rabbit mAb

Catalog No: #56868

Package Size: #56868-1 50ul #56868-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Complement factor B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Complement factor B Antibody detects endogenous levels of total Complement factor B
Immunogen Description	A synthesized peptide derived from human Complement factor B
Other Names	C3/C5 convertase; Complement factor B; PBF2; Properdin factor B;
Accession No.	Uniprot:P00751
Uniprot	P00751
Calculated MW	100kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of Complement factor B expression in Human plasma lysate.

Product Description

Factor B which is part of the alternate pathway of the complement system is cleaved by factor D into 2 fragments: Ba and Bb. Bb, a serine protease, then combines with complement factor 3b to generate the C3 or C5 convertase.

Note: This product is for in vitro research use only