

PROK1 Rabbit mAb

Catalog No: #56876

Package Size: #56876-1 50ul #56876-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PROK1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse
Specificity	PROK1 Antibody detects endogenous levels of total PROK1
Immunogen Description	A synthesized peptide derived from human PROK1
Other Names	EGVEGF; Mambakine; PK1; PRK1; PROK1; Prokineticin-1;
Accession No.	Uniprot:P58294
Uniprot	P58294
Calculated MW	14kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000FC:1:50

Images



Western blot analysis of PROK1 expression in Raw264.7 cell lysate.

Product Description

Potently contracts gastrointestinal (GI) smooth muscle. Induces proliferation, migration and fenestration (the formation of membrane discontinuities) in capillary endothelial cells derived from endocrine glands.

Note: This product is for in vitro research use only