

## TFF2 Rabbit mAb

Catalog No: #56880

Package Size: #56880-1 50ul #56880-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TFF2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	TFF2 Antibody detects endogenous levels of total TFF2
Immunogen Description	A synthesized peptide derived from human TFF2
Other Names	SML1; SP; Spasmolysin; TFF2; Trefoil factor 2;
Accession No.	Uniprot:Q03403
Uniprot	Q03403
Calculated MW	16kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

## Application Details

WB: 1:500~1:2000

## Images



Western blot analysis of TFF2 expression in Human stomach lysate.

## Product Description

Inhibits gastrointestinal motility and gastric acid secretion. Could function as a structural component of gastric mucus, possibly by stabilizing glycoproteins in the mucus gel through interactions with carbohydrate side chains.

---

Note: This product is for in vitro research use only