

mGluR3 Rabbit mAb

Catalog No: #56883

Package Size: #56883-1 50ul #56883-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	mGluR3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	mGluR3 Antibody detects endogenous levels of total mGluR3
Immunogen Description	A synthesized peptide derived from human mGluR3
Other Names	GLUR3; GPRC1C; GRM3; MGLu3; mGluR3;
Accession No.	Uniprot:Q14832
Uniprot	Q14832
Calculated MW	330kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500~1:2000

Images



Western blot analysis of mGluR3 expression in Human cerebellum lysate.

Product Description

Receptor for glutamate. The activity of this receptor is mediated by a G-protein that inhibits adenylate cyclase activity.

Note: This product is for in vitro research use only