

PCBP1 Rabbit mAb

Catalog No: #56893

Package Size: #56893-1 50ul #56893-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PCBP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	PCBP1 Antibody detects endogenous levels of total PCBP1
Immunogen Description	A synthesized peptide derived from human PCBP1
Other Names	Alpha-CP1; hnRNP E1; HNRPE1; HNRPX; PCBP1;
Accession No.	Uniprot:Q15365
Uniprot	Q15365
Calculated MW	40kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50

Images



Western blot analysis of PCBP1 expression in HeLa cell lysate.

Product Description

Single-stranded nucleic acid binding protein that binds preferentially to oligo dC.

Note: This product is for in vitro research use only