

CTNNA3 Rabbit mAb

Catalog No: #56900

Package Size: #56900-1 50ul #56900-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CTNNA3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	CTNNA3 Antibody detects endogenous levels of total CTNNA3
Immunogen Description	A synthesized peptide derived from human CTNNA3
Other Names	Alpha T catenin; Cadherin associated protein; Catenin alpha 3; Catna3; CTNNA3; VR22;
Accession No.	Uniprot:Q9UI47
Uniprot	Q9UI47
Calculated MW	100kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of CTNNA3 expression in 293 cell lysate.

Product Description

May be involved in formation of stretch-resistant cell-cell adhesion complexes.

Note: This product is for in vitro research use only