

p66 alpha Rabbit mAb

Catalog No: #56903

Package Size: #56903-1 50ul #56903-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	p66 alpha Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse
Specificity	p66 alpha Antibody detects endogenous levels of total p66 alpha
Immunogen Description	A synthesized peptide derived from human p66 alpha
Other Names	GATAD2A; Hp66alpha; P66 alpha; p66alpha;
Accession No.	Uniprot:Q86YP4
Uniprot	Q86YP4
Calculated MW	65,68kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200FC:1:50

Images



Western blot analysis of p66 alpha expression in K562 cell lysate.

Product Description

Transcriptional repressor.

Note: This product is for in vitro research use only