

EXOSC7 Rabbit mAb

Catalog No: #56917

Package Size: #56917-1 50ul #56917-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	EXOSC7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	EXOSC7 Antibody detects endogenous levels of total EXOSC7
Immunogen Description	A synthesized peptide derived from EXOSC7
Other Names	EAP1; Exosc7; p8; RRP42; Rrp42p;
Accession No.	Uniprot:Q15024
Uniprot	Q15024
Calculated MW	32kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of EXOSC7 expression in Jurkat cell lysate.

Product Description

Non-catalytic component of the RNA exosome complex which has 3'->5' exoribonuclease activity and participates in a multitude of cellular RNA processing and degradation events.

Note: This product is for in vitro research use only