

## Themis Rabbit mAb

Catalog No: #56918

Package Size: #56918-1 50ul #56918-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Themis Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	Themis Antibody detects endogenous levels of total Themis
Immunogen Description	A synthesized peptide derived from Themis
Other Names	Protein THEMIS; SPOT; THEMIS; TSEPA;
Accession No.	Uniprot:Q8N1K5
Uniprot	Q8N1K5
Calculated MW	73kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:50FC:1:50

## Images



Western blot analysis of Themis expression in Jurkat cell lysate.

## Product Description

Plays a central role in late thymocyte development by controlling both positive and negative T-cell selection. Required to sustain and/or integrate signals required for proper lineage commitment and maturation of T-cells.

---

Note: This product is for in vitro research use only