

## UCHL3 Rabbit mAb

Catalog No: #56922

Package Size: #56922-1 50ul #56922-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Product Name          | UCHL3 Rabbit mAb                                                                                   |
| Host Species          | Rabbit                                                                                             |
| Clonality             | Monoclonal                                                                                         |
| Isotype               | Rabbit IgG                                                                                         |
| Purification          | Affinity-chromatography                                                                            |
| Applications          | WB IHC                                                                                             |
| Species Reactivity    | Human Mouse Rat                                                                                    |
| Specificity           | UCHL3 Antibody detects endogenous levels of total UCHL3                                            |
| Immunogen Description | A synthesized peptide derived from UCHL3                                                           |
| Other Names           | Ubiquitin thioesterase L3; UCHL3;                                                                  |
| Accession No.         | Uniprot:P15374                                                                                     |
| Uniprot               | P15374                                                                                             |
| Calculated MW         | 26kDa                                                                                              |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle. |

## Application Details

WB:1:500~1:2000IHC:1:50~1:200

## Images



Western blot analysis of UCHL3 expression in K562 cell lysate.

## Product Description

Deubiquitinating enzyme (DUB) that controls levels of cellular ubiquitin through processing of ubiquitin precursors and ubiquitinated proteins. Thiol protease that recognizes and hydrolyzes a peptide bond at the C-terminal glycine of either ubiquitin or NEDD8.

---

Note: This product is for in vitro research use only