

UBE2B Rabbit mAb

Catalog No: #56925

Package Size: #56925-1 50ul #56925-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	UBE2B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	UBE2B Antibody detects endogenous levels of total UBE2B
Immunogen Description	A synthesized peptide derived from UBE2B
Other Names	hHR6B; HR6B; RAD6B; UBC2; UBE2B;
Accession No.	Uniprot:P63146
Uniprot	P63146
Calculated MW	17kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of UBE2B expression in Jurkat cell lysate.

Product Description

Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In association with the E3 enzyme BRE1 (RNF20 and/or RNF40), it plays a role in transcription regulation by catalyzing the monoubiquitination of histone H2B at 'Lys-120' to form H2BK120ub1.

Note: This product is for in vitro research use only