

GTF2F2 Rabbit mAb

Catalog No: #56926

Package Size: #56926-1 50ul #56926-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GTF2F2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	GTF2F2 Antibody detects endogenous levels of total GTF2F2
Immunogen Description	A synthesized peptide derived from GTF2F2
Other Names	BTF4; Gtf2f2; RAP30; TF2F2; TFIIF; TFIIF beta;
Accession No.	Uniprot:P13984
Uniprot	P13984
Calculated MW	30kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of GTF2F2 expression in HeLa cell lysate.

Product Description

TFIIF is a general transcription initiation factor that binds to RNA polymerase II and helps to recruit it to the initiation complex in collaboration with TFIIB. It promotes transcription elongation. This subunit shows ATP-dependent DNA-helicase activity.

Note: This product is for in vitro research use only