

AP2S1 Rabbit mAb

Catalog No: #56927

Package Size: #56927-1 50ul #56927-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	AP2S1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human Mouse Rat
Specificity	AP2S1 Antibody detects endogenous levels of total AP2S1
Immunogen Description	A synthesized peptide derived from AP2S1
Other Names	AP17; Ap2s1; CLAPS2; Sigma adaptin 3b; Sigma2 adaptin;
Accession No.	Uniprot:P53680
Uniprot	P53680
Calculated MW	17kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of AP2S1 expression in (1) 293 cell lysate; (2) Mouse brain lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Product Description

Component of the adaptor protein complex 2 (AP-2). Adaptor protein complexes function in protein Transport via Transport vesicles in different membrane traffic pathways.

Note: This product is for in vitro research use only