

## GBP1 Rabbit mAb

Catalog No: #56930

Package Size: #56930-1 50ul #56930-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	GBP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	GBP1 Antibody detects endogenous levels of total GBP1
Immunogen Description	A synthesized peptide derived from human GBP1
Other Names	GBP1; GTP binding protein 1; Guanine nucleotide binding protein 1; Guanylate binding protein 1; HuGBP1; Interferon induced guanylate binding protein 1;
Accession No.	Uniprot:P32455
Uniprot	P32455
Calculated MW	68kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

## Images



Western blot analysis of GBP1 expression in HeLa cell lysate.

## Product Description

Binds GTP, GDP and GMP.

---

Note: This product is for in vitro research use only