

IL8 Rabbit mAb

Catalog No: #56932

Package Size: #56932-1 50ul #56932-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IL8 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	IL8 Antibody detects endogenous levels of total IL8
Immunogen Description	A synthesized peptide derived from human IL8
Other Names	CXCL8; GCP1; Granulocyte chemotactic protein 1; IL8; Inteleukin 8; LECT; LUCT; LY NAP; MDNCF; MONAP; NAF; NAP1; Neutrophil activating protein 1; SCYB8; TSG1;
Accession No.	Uniprot:P10145
Uniprot	P10145
Calculated MW	11kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of IL8 expression in IL8 recombinant protein cell lysate.

Product Description

IL-8 is a chemotactic factor that attracts neutrophils, basophils, and T-cells, but not monocytes. It is also involved in neutrophil activation.

Note: This product is for in vitro research use only