

NUDT5 Rabbit mAb

Catalog No: #56936

Package Size: #56936-1 50ul #56936-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NUDT5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	NUDT5 Antibody detects endogenous levels of total NUDT5
Immunogen Description	A synthesized peptide derived from human NUDT5
Other Names	hYSAH1; Nudt5; YSA1; YSA1H;
Accession No.	Uniprot:Q9UUKK9
Uniprot	Q9UUKK9
Calculated MW	30kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of NUDT5 expression in HeLa cell lysate.

Product Description

Hydrolyzes with similar activities ADP-ribose, ADP-mannose, and ADP-glucose. Can also hydrolyze other nucleotide sugars with low activity.

Note: This product is for in vitro research use only